

# AdvanceBio HIC: a Hydrophobic HPLC Column for Monoclonal Antibody (mAb) Variant Analysis

Using the Agilent 1260 Infinity II Bio-Inert LC

#### **Authors**

Andrew Coffey and Sandeep Kondaveeti Agilent Technologies, Inc.

# **Abstract**

This Application Note describes the separation of oxidized monoclonal antibody (mAb) variants from their native form using the Agilent AdvanceBio HIC column. Oxidation of exposed amino acid side chain residues such as methionine, cysteine, and tryptophan is a common degradation pathway for monoclonal antibodies, and presents a major analytical challenge in biotechnology. Often, oxidized mAbs have decreased potency compared to their native form<sup>1</sup>. Therefore, to ensure the therapeutic efficacy of the mAb products, analysis of such degradation is critical. Oxidation of amino acid residues on an mAb can alter the hydrophobic nature of the mAb by the increase in polarity of the oxidized form, or also due to resulting conformational changes<sup>2</sup>. HPLC methods for separating biomolecules based on differences in hydrophobicity include reversed-phase and hydrophobic interaction chromatography (HIC). HIC can be applied to characterize mAb variants resulting from post-translational modifications (PTMs). The AdvanceBio HIC column provides excellent resolution of oxidized mAb variants from unmodified forms, and can resolve oxidized species without mAb digestion into subunits or other sample preparation methods.

# Introduction

mAbs and related products such as antibody drug conjugates (ADCs) and bispecific antibodies (bsAbs) are the fastest growing classes of biotherapeutics. Recombinant mAbs are subject to many PTMs during processing, delivery, and storage. Among these modifications, oxidation of exposed amino acid side chains such as methionine (Met) and tryptophan (Trp) is a common occurrence. Various researchers have reported that oxidation of mAbs has an adverse effect on product shelf life and bio-activity<sup>1,2</sup>. Therefore, developing analytical methods to detect oxidized mAb variants has gained interest. The sulfoxide and sulfone side chains of methionine-oxidized mAb products are larger and more polar compared to the native form, which may alter protein structure, stability, and biological function. Hydrophobicity-based HPLC methods, such as reversed-phase liquid chromatography (RPLC) and HIC, are often used to characterize mAb variants. Recently, several studies have indicated that HIC can be applied to monitor oxidation of recombinant mAbs with reasonable selectivity and ease, as an excellent alternative to RPLC3.

HIC is similar to RPLC in that separation of analytes is based on hydrophobic interactions with the stationary phase. The elution order in HIC enables proteins to be ranked based on their relative hydrophobicity. Unlike RPLC, HIC employs nondenaturing conditions, does not require the use of organic solvents or high temperatures, and separations

are carried out at physiological pH, allowing for the preservation of protein structure. Thus, conformational changes in the native form of the protein may be analyzed using HIC<sup>4</sup>.

AdvanceBio HIC is a silica-based HPLC column designed for the separation of mAbs and related products. Its unique proprietary bonded phase chemistry provides high resolution and desired selectivity for the analysis of mAbs and mAb variants. This Application Note describes the separation of oxidized NIST mAb variants using an AdvanceBio HIC column.

# **Experimental**

## **Equipment and materials**

All chemicals and reagents were HPLC grade or higher, and were obtained from Sigma-Aldrich (now Merck) or VWR Scientific. Humanized IgG1k mAb sample (product item no. 8671) was obtained from NIST SRM Standards. Water was purified using a Milli-Q A10 water purification system (Millipore).

### Instrumentation

Agilent 1260 Infinity II bio-inert LC comprising:

- Agilent 1260 Infinity II bio-inert pump (G5654A)
- Agilent 1260 Infinity II bio-inert multisampler (G5668A) with sample cooler (option no. 100)
- Agilent 1260 Infinity II multicolumn thermostat (G7116A) with bio-inert heat exchanger (option no. 019)
- Agilent 1260 Infinity II diode array detector WR (G7115A) with bio-inert flow cell (option no. 028)

#### Software

Agilent OpenLab 2.2 CDS

## mAb Oxidation with t-BHP treatment

A solution of 1 mL of NIST mAb (1 mg/mL) in formulation buffer was incubated with 0.2 % (v/v) of 70 % tert-butyl hydroperoxide (t-BHP) solution at room temperature for 24 hours. Residual oxidant was removed using Amicon Ultra-10 centrifugal filters with a molecular weight cut-off of 10 kDa. Samples were then buffer exchanged with 50 % mobile phase B.

Reaction conditions used to obtain Figure 6 data: 2 % (v/v) of 70 % *t*-BHP solution was added to a 1-mL sample of NIST mAb (1 mg/mL), and the reaction mixture was injected onto the column. The sample vial was held at 7 °C, and multiple injections from the same vial were carried out.

## mAb Oxidation with H<sub>2</sub>O<sub>2</sub> treatment

A solution of 1 mL of NIST mAb (1 mg/mL) in formulation buffer was incubated with 0.2 % (v/v) of 50 % hydrogen peroxide ( $H_2O_2$ ) solution at room temperature for 24 hours. Residual oxidant was removed using Amicon Ultra-10 centrifugal filters with a molecular weight cut-off of 10 kDa. Samples were then buffer exchanged with 50 % mobile phase B.

# **Results and discussion**

Protein oxidation is frequently monitored in stability studies or stressed samples during formulation development (for example, samples exposed to a chemical oxidant such as hydrogen peroxide (H2O2), UV light, or metal ions). In this study, t-BHP and H<sub>2</sub>O<sub>2</sub> were used as chemical oxidants to promote oxidation of NIST mAb samples. It was previously reported that both of these reagents tend to specifically oxidize Met side chain residues of the mAb. H<sub>2</sub>O<sub>2</sub> more readily oxidizes less accessible, buried residues. whereas t-BHP is known to target more surface-exposed Met residues7. Figure 1 illustrates the reaction scheme for Met oxidation induced by chemical oxidants.

**Figure 1.** Methionine oxidation induced by chemical oxidant.

The NIST mAb (humanized IgG1k) amino acid sequence in Figure 2 shows that there are six possible surface-accessible Met residues located on both heavy chains of the mAb. Based on prior studies for most human IgG1-subclass antibodies, Met residues localized to the CH<sub>2</sub> and CH<sub>3</sub> domains of the antigen binding, or Fc, region are known to be highly susceptible to oxidation<sup>5</sup>. In the case of NIST mAb, Met 255 and Met 431 correspond to the amino acid residues prone to oxidation. This is depicted by the illustration in Figure 3.

## Method conditions

	HPLC conditions			
Column	AdvanceBio HIC, 4.6 × 100 mm (p/n 685975-908)			
Mobile phase	Eluent A) 50 mM sodium phosphate, pH 7.0 Eluent B) 2 M ammonium sulfate in 50 mM sodium phosphate, pH 7.0			
Flow rate	0.3 to 0.5 mL/min			
Column temperature	25 °C			
Injection volume	5 μL			
Final sample concentration	1 mg/mL			
Detection	UV, 220 nm			
Gradient profile	Flow rate: 0.5 mL/min  Time %A %B 0 50 50 20 100 0 25 100 0 30 50 50 40 50 50			
Shallower gradient	Flow rate: 0.3 mL/min  Time %A %B 0 40 60 40 90 10 45 90 10 50 40 60 60 40 60			

#### Heavy chain

QVTLRESGPA	LVKPTQTLTL	TCTFSGFSLS	TAG <b>M</b> SVGWIR	QPPGKALEWL	ADIWWDDKKH	YNPSLKDRLT	
ISKDTSKNQV	${\tt VLKVTN}{\bf M}{\tt DPA}$	DTATYYCARD	MIFNFYFDVW	GQGTTVTVSS	ASTKGPSVFP	LAPSSKSTSG	
GTAALGCLVK	DYFPEPVTVS	WNSGALTSGV	HTFPAVLQSS	GLYSLSSVVT	VPSSSLGTQT	YICNVNHKPS	
NTKVDKRVEP	KSCDKTHTCP	PCPAPELLGG	PSVFLFPPKP	KDTL <b>m</b> isrtp	EVTCVVVDVS	HEDPEVKFNW	
YVDGVEVHNA	KTKPREEQYN	STYRVVSVLT	VLHQDWLNGK	EYKCKVSNKA	LPAPIEKTIS	KAKGQPREPQ	
VYTLPPSREE	<b>M</b> TKNQVSLTC	LVKGFYPSDI	AVEWESNGQP	ENNYKTTPPV	LDSDGSFFLY	SKLTVDKSRW	
OOGNVFSCSV	MHEALHNHYT	OKSLSLSPGK					

### Light chain

DIQMTQSPST	LSASVGDRVT	ITCSASSRVG	Y <b>M</b> HWYQQKPG	KAPKLLIYDT	SKLASGVPSR	FSGSGSGTEF
TLTISSLQPD	DFATYYC <mark>FQG</mark>	SGYPFTFGGG	TKVEIKRTVA	APSVFIFPPS	DEQLKSGTAS	VVCLLNNFYP
REAKVQWKVD	NALQSGNSQE	SVTEQDSKDS	TYSLSSTLTL	SKADYEKHKV	YACEVTHQGL	SSPVTKSFNR GEC

Figure 2. NIST mAb amino acid sequence.




Figure 3. Methionine residues located in Fc region are most susceptible to oxidation in human IgG1 mAbs.

An AdvanceBio HIC column was able to differentiate oxidized mAb variants from the untreated mAb sample under low salt starting conditions. Oxidation of the NIST mAb with t-BHP under reported experimental conditions resulted in multiple peaks with shorter retention times, presumably due to conformational change. The HIC chromatogram (Figure 4) showing earlier retained peaks labeled 1 to 6 likely indicates the result of oxidized Met residues on the mAb, and peak 7 with a retention time of approximately 12.6 minutes, corresponds to nonoxidized mAb. For the H<sub>2</sub>O<sub>2</sub>-treated mAb sample, complete oxidation occurred, with three peaks eluting in a shorter retention time, indicating more aggressive oxidation of Met residues. These differences in the chromatograms of the IgG1k mAb sample incubated with two different oxidation reagents suggest that reactivity is governed by solvent accessibility of the Met residues and steric limitations of the oxidizing agent, as previously reported6.

To further improve the resolution, a slower and shallower gradient was used. Using a flow rate of 0.3 mL/min and a starting ammonium sulfate concentration of 1.2 M with a lower gradient rate of 25 mM/min, better resolution was achieved with a relatively short analysis time (Figure 5). In this chromatogram, multiple mAb-oxidized species are clearly observed from the untreated mAb sample.




Figure 4. Separation of oxidized NIST mAb variants using lower starting salt concentration.





Figure 5. Separation of oxidized NIST mAb variants using a shallow gradient.

In Figure 6, the NIST mAb sample was incubated with 2 % (v/v) t-BHP, and the oxidation reaction was monitored at various time points using shallower gradient conditions. As represented by an overlay of chromatograms, the mAb oxidation progressed with t-BHP incubation time. Multiple mAb oxidation species were observed within a few hours of the oxidation reaction. This suggested that surface-accessible Met residues in both heavy chains of the mAb sample might be oxidized randomly,

which was previously reported<sup>5</sup>. Further oxidation of the mAb sample after 10 hours of reaction led to a broad peak, indicating forced oxidation. It has previously been speculated that oxidation of deeply buried Met residues can lead to a more dramatic structural change, which may cause the mAb to partially unfold<sup>7</sup>. Partially unfolded mAb is likely to have more conformational variation, resulting in a broader peak with a large retention time shift.

# **Conclusions**

The AdvanceBio HIC column demonstrated the separation of oxidized mAb variants from its native form. Using the AdvanceBio HIC column, optimal separation of oxidized mAb variants can be achieved using slower flow rates and shallower gradient conditions, while maintaining relatively short analysis times



**Figure 6.** Monitoring the *t*-BHP oxidized mAb reaction.

Time	%A	%B	
0	40	60	25 mM/min
40	90	10	23 111101/111111
45	90	10	
50	40	60	
60	40	60	

# References

- 1. Zhang, Y.; et al. Hydrophobic interaction chromatography of soluble Interleukin I receptor type II to reveal chemical degradations resulting in loss of potency, Anal. Chem. **2008**, 80(18), 7022-8.
- 2. Gaza-Bulseco, G.; et al. Effect of methionine oxidation of a recombinant monoclonal antibody on the binding affinity to protein A and protein G. J. Chromatogr. B Analyt. Technol. Biomed. Life Sci. **2008**, 870(1), 56–62.
- 3. Boyd, D.; Kaschak, T.; Yan, B. HIC resolution of an IgG1 with an oxidized Trp in a complementarity determining region, *J. Chromatogr. B.* **2011**, *879*(*13*–*14*), 955-60.
- Fekete, S.; et al. Hydrophobic interaction chromatography for the characterization of monoclonal antibodies and related products, J. Pharm. Biomed. Anal. 2016, 130, 3–18.

- Chumsae, C.; et al. Comparison of methionine oxidation in thermal stability and chemically stressed samples of a fully human monoclonal antibody. J. Chromatogr. B: Analyt. Technol. Biomed. Life Sci. 2007, 850(1-2), 285-294.
- Shen, F. J.; et al. The application of tert-butylhydroperoxide oxidation to study sites of potential methionine oxidation in a recombinant antibody. In *Techniques in* Protein Chemistry; Marshak, D. R., Ed.; Academic Press, Inc.: San Diego, CA, **1996**; 7, 275–284.
- 7. Liu, D.; et al. Structure and stability changes of human IgG1 Fc as a consequence of methionine oxidation. *Biochemistry* **2008**, *47*(*18*), 5088–5100.

## www.agilent.com/chem

This information is subject to change without notice.

For Research Use Only. Not for use in diagnostic procedures.

