

Interpretation Manual - Gastric or Gastroesophageal Junction Adenocarcinoma

PD-L1 IHC 22C3 pharmDx is FDA-approved for in vitro diagnostic use




Table of Contents

Intended Use	04
Introduction	06
PD-L1 Overview	08
PD-L1 IHC 22C3 pharmDx Overview	10
Kit Configuration (SK006)	11
Technical Considerations	12
Specimen Preparation	12
In-house Control Tissue	12
Optional Additional In-house Control: Tonsil Tissue	13
Tissue Processing	13
PD-L1 IHC 22C3 pharmDx Staining Procedure	14
Technical Checklist	17
Slide Evaluation	18
General Considerations	18
Specimen Adequacy	18
Evaluating Controls	18
Slide Evaluation Flow Chart	22
Combined Positive Score	23
Definition of Combined Positive Score (CPS)	23
CPS Numerator Inclusion and Exclusion Criteria	23
Determining Combined Positive Score	23
Suggested Methods	25
Interpretation of CPS	28
Identifying Patients with Gastric or GEJ Adenocarcinoma for Treatment	29
PD-L1 IHC 22C3 pharmDx Testing Scheme	30
Reporting Results	31
Combined Positive Score Summary and Examples	32
Key Considerations in Scoring PD-L1 IHC 22C3 pharmDx Stained Specimens	32
Image Guide for Interpretation of PD-L1 IHC 22C3 pharmDx Staining in Gastric or GEJ Adenocarcinoma	33
CPS < 1 Case Examples	47
CPS ≥ 1 Case Examples	53
Near Cut-off Case Examples (CPS 0–10)	59
Artifacts	62
Troubleshooting Guide	67
Clinical Performance Evaluation	68
References	70

Intended Use

For in vitro diagnostic use.

PD-L1 IHC 22C3 pharmDx is a qualitative immunohistochemical assay using Monoclonal Mouse Anti-PD-L1, Clone 22C3 intended for use in the detection of PD-L1 protein in formalin-fixed, paraffin-embedded (FFPE) non-small cell lung cancer (NSCLC) and gastric or gastroesophageal junction (GEJ) adenocarcinoma tissues using EnVision FLEX visualization system on Autostainer Link 48.

Non-Small Cell Lung Cancer (NSCLC)


PD-L1 protein expression in NSCLC is determined by using Tumor Proportion Score (TPS), which is the percentage of viable tumor cells showing partial or complete membrane staining at any intensity. The specimen should be considered to have PD-L1 expression if $TPS \geq 1\%$ and high PD-L1 expression if $TPS \geq 50\%$.

PD-L1 IHC 22C3 pharmDx is indicated as an aid in identifying NSCLC patients for treatment with KEYTRUDA® (pembrolizumab). See the KEYTRUDA product label for expression cutoff values guiding therapy in specific clinical circumstances.

Gastric or Gastroesophageal Junction (GEJ) Adenocarcinoma

PD-L1 protein expression in gastric or GEJ adenocarcinoma is determined by using Combined Positive Score (CPS), which is the number of PD-L1 staining cells (tumor cells, lymphocytes, macrophages) divided by the total number of viable tumor cells, multiplied by 100. The specimen should be considered to have PD-L1 expression if $CPS \geq 1$.

PD-L1 IHC 22C3 pharmDx is indicated as an aid in identifying gastric or GEJ adenocarcinoma patients for treatment with KEYTRUDA (pembrolizumab).



Introduction

PD-L1 IHC 22C3 pharmDx is the only companion diagnostic approved by the FDA as an aid in identifying patients with gastric or gastroesophageal junction (GEJ) adenocarcinoma for treatment with KEYTRUDA® (pembrolizumab). This Interpretation Manual is provided as a tool to help guide pathologists and laboratory personnel in achieving correct and reproducible results in assessing PD-L1 expression in formalin-fixed, paraffin-embedded gastric or GEJ adenocarcinoma specimens. PD-L1 expression evaluation may be used to identify patients for anti-PD-1 immunotherapy.

The manual provides detailed scoring guidelines and technical information from the PD-L1 IHC 22C3 pharmDx Instructions for Use (IFU) to ensure high-quality staining and diagnostic assessment. To help familiarize you with the requirements for scoring gastric or GEJ adenocarcinoma stains with PD-L1 IHC 22C3 pharmDx, example cases of various PD-L1 expression levels are provided as reference. These example cases and in-depth recommendations for interpretation of gastric or GEJ adenocarcinoma specimens stained with PD-L1 IHC 22C3 pharmDx can help individual labs achieve reproducible and reliable results.

PD-L1 IHC 22C3 pharmDx is considered a qualitative immunohistochemistry assay. PD-L1 expression in gastric or GEJ adenocarcinoma is determined by using Combined Positive Score (CPS), which is the number of PD-L1 staining cells (tumor cells, lymphocytes, macrophages) divided by the total number of viable tumor cells, multiplied by 100.

Gastric or GEJ adenocarcinoma tissue specimens that are tested for PD-L1 expression are scored and divided into two groups based on their Combined Positive Score (CPS):

- CPS < 1: No PD-L1 expression
- CPS ≥ 1: PD-L1 expression

For more details on staining and interpretation, please refer to the current version of the IFU provided with PD-L1 IHC 22C3 pharmDx, Code SK006 or visit www.agilent.com.

Assay Interpretation

The clinical interpretation of any staining, or the absence of staining, must be complemented by the evaluation of proper controls. Evaluation must be made by a qualified pathologist within the context of the patient's clinical history and other diagnostic tests. This product is intended for in vitro diagnostic (IVD) use.

Reporting Results

To help understand what information should be reported to the treating physician, please refer to the Reporting Results section of this manual on page 26.

Photomicrographs

The included photomicrographs are of gastric or GEJ adenocarcinoma unless otherwise noted.

Note: Photomicrograph magnification levels may appear different than indicated in respective annotations due to adjustment of image size.

PD-L1 Overview

The PD-1/PD-L1 Pathway Controls the Immune Response in Normal Tissue

Programmed death-ligand 1 (PD-L1) is a transmembrane protein that binds to the programmed death-1 receptor (PD-1) during immune system modulation. The PD-1 receptor is typically expressed on cytotoxic T-cells and other immune cells, while the PD-L1 ligand is typically expressed on normal cells. Normal cells use the PD-1/PD-L1 interaction as a mechanism of protection against immune recognition by inhibiting the action of T-cells (Figure 1). Inactivation of cytotoxic T-cells downregulates the immune response such that the inactive T-cell is exhausted, ceases to divide, and might eventually die by programmed cell death, or apoptosis.

The Tumor Escapes Detection by Utilizing the PD-1/PD-L1 Pathway

Many tumor cells are able to upregulate the expression of PD-L1 as a mechanism to evade the body's natural immune response. Activated T-cells recognize the PD-L1 marker on the tumor cell, similar to that of a normal cell, and PD-L1 signaling renders the T-cell inactive (Figure 2). The tumor cell escapes the immune cycle, continues to avoid detection for elimination, and is able to proliferate.

Anti-PD-1 Therapy Enables the Immune Response Against Tumors

PD-1/PD-L1 interaction between tumor cells and activated T-cells (Figure 3) is a mechanistic pathway used by immunotherapeutic agents. When the tumor cell is unable to interact with the activated T-cell, the immune system remains active, helping to prevent immunosuppression.

PD-L1 IHC 22C3 pharmDx Detects PD-L1 in Urothelial Carcinoma Specimens

Detection of PD-L1 upregulation in urothelial carcinoma is a biomarker for response to anti-PD-1 therapy. PD-L1 IHC 22C3 pharmDx is the only companion diagnostic used in the KEYTRUDA® (pembrolizumab) clinical trial (KEYNOTE-052) to evaluate the relationship between PD-L1 expression and clinical efficacy. KEYTRUDA is a humanized monoclonal PD-1-blocking antibody.

Note: Gastric cancer specimens include gastric or gastroesophageal junction adenocarcinoma.




Figure 1: Inactivation of T-cells limits damage to normal tissue.




Figure 2: Inactivation of T-cells reduces tumor cell death and elimination.




Figure 3: Blocking the PD-1/PD-L1 interaction helps to enable active T-cells and tumor cell death and elimination.

PD-L1 IHC 22C3 pharmDx Overview

What is PD-L1 IHC 22C3 pharmDx?

PD-L1 IHC 22C3 pharmDx is the only companion diagnostic indicated as an aid in identifying patients with gastric or gastroesophageal junction (GEJ) adenocarcinoma for treatment with KEYTRUDA® (pembrolizumab). PD-L1 IHC 22C3 pharmDx is a qualitative immunohistochemical (IHC) assay intended for use in the detection of PD-L1 protein in formalin-fixed, paraffin-embedded (FFPE) gastric cancer tissue samples using Autostainer Link 48.

Components of PD-L1 IHC 22C3 pharmDx

PD-L1 IHC 22C3 pharmDx contains optimized reagents to perform an IHC staining procedure using a linker and a chromogen enhancement reagent (Figure 4). Deparaffinization, rehydration, and target retrieval is performed using a 3-in-1 procedure on PT Link. Following peroxidase block, specimens are incubated with monoclonal mouse primary antibody to PD-L1 or the Negative Control Reagent. Specimens are then incubated with a Mouse LINKER, followed by incubation with a ready-to-use Visualization Reagent consisting of secondary antibody molecules and horseradish peroxidase molecules coupled to a dextran polymer backbone.

The enzymatic conversion of the subsequently added chromogen results in precipitation of a visible reaction product at the site of the antigen. The color of the chromogenic reaction is modified by a chromogen enhancement reagent. The specimen may then be counterstained and coverslipped. Results are interpreted using a light microscope.




Figure 4: PD-L1 IHC 22C3 pharmDx staining procedure.

Kit Configuration




Figure 5: PD-L1 IHC 22C3 pharmDx components.

PD-L1 IHC 22C3 pharmDx (Code SK006) contains reagents to perform 50 tests in up to 15 individual runs (Figure 5):

- 1 EnVision FLEX Target Retrieval Solution, Low pH (50×)
- 2 Peroxidase-Blocking Reagent
- 3 Primary antibody: Monoclonal Mouse Anti-PD-L1, Clone 22C3
- 4 Negative Control Reagent
- 5 Mouse LINKER
- 6 Visualization Reagent-HRP
- 7 DAB+ Substrate Buffer
- 8 DAB+ Chromogen
- 9 DAB Enhancer
- 10 PD-L1 IHC 22C3 pharmDx Control Cell Line Slides*

* Dr. AF Gazdar and Dr. JD Minna at NIH are acknowledged for their contribution in developing NCI-H226 (ATCC Number: CRL-5826™)

EnVision FLEX Wash Buffer, (20×) (Code K8007) and EnVision FLEX Hematoxylin (Code K8008) are required but not included in the kit.



Technical Considerations

Technical problems related to PD-L1 IHC 22C3 pharmDx may arise and can be attributed to two factors: specimen collection and preparation prior to performing the test, and the actual performance of the test itself. Technical problems are generally related to procedural deviations and can be controlled and minimized through training and, where necessary, clarification of the product instructions.

Specimen Preparation

Specimens must be handled to preserve the tissue for immunohistochemical staining. Determine intact tumor morphology and the presence of sufficient tumor cells for evaluation. Use standard methods of tissue processing for all specimens.

In-house Control Tissue

Differences in processing and embedding in the user's laboratory may produce significant variability in results. Include positive and negative in-house control tissue in each staining run, in addition to the PD-L1 IHC 22C3 pharmDx Control Cell Line Slide.

Select positive and negative control tissue from fresh specimens of the same tumor indication as the patient specimen. Fix, process, and embed the control tissue in the same manner. Control tissues processed differently from the patient specimen validate reagent performance only and do not verify tissue preparation.

The ideal positive control tissue provides a complete dynamic representation of weak to moderate staining of tumor cells and tumor-associated mononuclear inflammatory cells (MICs: lymphocytes and macrophages). The negative control tissue should demonstrate a CPS less than 1 ($CPS < 1$); with no staining in tumor cells and only a few staining immune cells.

Optional Additional In-house Control: Tonsil Tissue

Tonsil stained with PD-L1 should be pre-screened to exhibit strong staining in portions of the crypt epithelium and weak to moderate staining of the follicular macrophages in the germinal centers. PD-L1 expression of the endothelium, fibroblasts, as well as the surface epithelium should be negative.

Tissue Processing

Formalin-fixed, paraffin-embedded tissues have been validated for use. Block specimens into a thickness of 3 mm or 4 mm, fix in formalin and dehydrate and clear in a series of alcohols and xylene, followed by infiltration with melted paraffin. The paraffin temperature should not exceed 60 °C. Feasibility studies on NSCLC tissue samples were performed with fixation in 10% neutral buffered formalin for 12–72 hours. Fixation times of 3 hours or less should not be used for PD-L1 assessment. The use of PD-L1 IHC 22C3 pharmDx on decalcified tissues or tissues processed with other fixatives has not been validated and is not recommended.

Cut tissue specimens into sections of 4–5 µm. After sectioning, tissues should be mounted on Dako FLEX IHC microscope slides (Code K8020) or Fisherbrand Superfrost Plus slides, and then placed in a 58 ± 2 °C oven for 1 hour. Store tissue sections in the dark at 2–8 °C (preferred) or at room temperature in the dark up to 25 °C to preserve antigenicity, and stain within 6 months of sectioning.

PD-L1 IHC 22C3 pharmDx Staining Procedure

The PD-L1 IHC 22C3 pharmDx reagents and instructions have been designed for optimal performance. Further dilution of the reagents, alteration of incubation times, temperatures, or materials may give erroneous results. All of the required steps and incubation times for staining are pre-programmed in the DakoLink software.

Reagent Storage

Store all components of PD-L1 IHC 22C3 pharmDx, including Control Cell Line Slides, in the dark at 2–8 °C when not in use.

Reagent Preparation

Equilibrate all components to room temperature (20–25 °C) prior to immunostaining. Do not use after the expiration date printed on the outside of the package.

EnVision FLEX Target Retrieval Solution, Low pH

Dilute EnVision FLEX Target Retrieval Solution, Low pH, (50x) 1:50 using distilled or deionized water (reagent-quality water). One 30 mL bottle of concentrate provides 1.5 L of working solution, which is sufficient to fill one PT Link tank. Discard 1x EnVision FLEX Target Retrieval Solution, Low pH after 3 uses or 5 days after dilution.

EnVision FLEX Wash Buffer

Dilute EnVision FLEX Wash Buffer, (20x) 1:20 using distilled or deionized water (reagent-quality water). Store unused 1x buffer at 2–8 °C for no more than one month. Discard if cloudy in appearance.

DAB+ Substrate-Chromogen Solution

Add 1 drop of DAB+ Chromogen per mL of DAB+ Substrate Buffer and mix. Prepared DAB+ Substrate-Chromogen is stable for 5 days if stored in the dark at 2–8 °C. Mix the DAB+ Substrate-Chromogen Solution thoroughly prior to use. Any precipitate developing in the solution will not affect staining quality.

- *If using an entire bottle of DAB+ Substrate Buffer, add 9 drops of DAB+ Chromogen. Although the DAB+ Substrate Buffer label states 7.2 mL, this is the usable volume and does not account for the “dead volume” of DAB+ Substrate Buffer in the bottle.*
- The color of the DAB+ Chromogen may vary from clear to lavender brown. This will not affect the performance of the product. Dilute per the guidelines above. Adding excess DAB+ Chromogen to the DAB+ Substrate Buffer results in deterioration of the positive signal

Controls to Assess Staining Quality

The following quality controls should be included in each staining run:

- One PD-L1 IHC 22C3 pharmDx Control Cell Line Slide stained with the primary antibody
- Positive and negative in-house control tissues stained with the primary antibody
- Subsequent sections of each patient specimen stained with the Negative Control Reagent

Deparaffinization, Rehydration, and Target Retrieval

Use PT Link to perform a Deparaffinization, Rehydration, and Target Retrieval 3-in-1 procedure.

- Set Preheat and Cool to 65 °C, and set Heat to 97 °C for 20 minutes
- Fill PT Link tanks with 1.5 L per tank of 1x EnVision FLEX Target Retrieval Solution, Low pH, working solution to cover the tissue sections
- Preheat the Target Retrieval Solution, Low pH to 65 °C
- Immerse Autostainer racks containing mounted, FFPE tissue sections into the preheated Target Retrieval Solution, Low pH in PT Link tank. Incubate for 20 minutes at 97 °C
- When incubation has been completed and the temperature has cooled to 65 °C, remove each Autostainer slide rack with slides from the PT Link tank and immediately place the slides into a tank (e.g., PT Link Rinse Station, Code PT109) containing room temperature 1x EnVision FLEX Wash Buffer working solution
- Leave Autostainer rack with slides in room temperature 1x EnVision FLEX Wash Buffer for 5 minutes

Staining and Counterstaining

- Place the Autostainer rack with slides on the Autostainer Link 48
- Ensure slides remain wet with buffer while loading and prior to initiating the run. Dried tissue sections may display increased non-specific staining
- Select the PD-L1 IHC 22C3 pharmDx protocol. The instrument performs the staining and counterstaining procedures by applying the appropriate reagent, monitoring the incubation time, and rinsing slides between reagents
- Counterstain slides using EnVision FLEX Hematoxylin, Code K8008

Mounting

Use non-aqueous permanent mounting media. To minimize fading, store slides in the dark at room temperature (20–25 °C).

Technical Checklist

Use the checklist below to ensure correct usage of PD-L1 IHC 22C3 pharmDx:

Customer Name/Institution _____

Name and Title _____

Autostainer Link 48 Serial Number _____ Software Version _____

	Yes	No
Regular preventive maintenance is performed on the Autostainer Link 48 and PT Link?	<input type="checkbox"/>	<input type="checkbox"/>
PD-L1 IHC 22C3 pharmDx is used before the expiration date printed on the outside of the box?	<input type="checkbox"/>	<input type="checkbox"/>
All PD-L1 IHC 22C3 pharmDx components, including Control Cell Line Slides, are stored in the dark at 2–8 °C?	<input type="checkbox"/>	<input type="checkbox"/>
All PD-L1 IHC 22C3 pharmDx components, including Control Cell Line Slides, are equilibrated to room temperature (20–25 °C) prior to immunostaining?	<input type="checkbox"/>	<input type="checkbox"/>
Appropriate positive and negative control tissue from urothelial carcinoma are identified?	<input type="checkbox"/>	<input type="checkbox"/>
Tissues are fixed in neutral buffered formalin?	<input type="checkbox"/>	<input type="checkbox"/>
Tissues are infiltrated with melted paraffin, at or below 60 °C?	<input type="checkbox"/>	<input type="checkbox"/>
Tissue sections of 4–5 µm are mounted on Dako FLEX IHC Microscope Slides or Fisherbrand Superfrost Plus charged slides?	<input type="checkbox"/>	<input type="checkbox"/>
Specimens are oven-dried at 58 ± 2 °C for 1 hour?	<input type="checkbox"/>	<input type="checkbox"/>
Specimens are stained within 1 month of sectioning when stored in the dark at 2–8 °C (preferred) or at room temperature in the dark up to 25 °C?	<input type="checkbox"/>	<input type="checkbox"/>
EnVision FLEX Target Retrieval Solution, Low pH is prepared properly? pH of 1× Target Retrieval Solution must be 6.1 ± 0.2.	<input type="checkbox"/>	<input type="checkbox"/>
EnVision FLEX Wash Buffer is prepared properly?	<input type="checkbox"/>	<input type="checkbox"/>
DAB+ Substrate-Chromogen Solution is prepared properly?	<input type="checkbox"/>	<input type="checkbox"/>
Slides are counterstained with EnVision FLEX Hematoxylin?	<input type="checkbox"/>	<input type="checkbox"/>
The Deparaffinization, Rehydration, and Target Retrieval 3-in-1 procedure is followed using PT Link?	<input type="checkbox"/>	<input type="checkbox"/>
Slides remain wet with buffer while loading and prior to initiating run on Autostainer Link 48?	<input type="checkbox"/>	<input type="checkbox"/>
The PD-L1 IHC 22C3 pharmDx protocol is selected on Autostainer Link 48?	<input type="checkbox"/>	<input type="checkbox"/>
Do you have all the necessary equipment to perform the PD-L1 IHC 22C3 pharmDx according to protocol? If not, specify what is missing in comments below.	<input type="checkbox"/>	<input type="checkbox"/>

Additional observations or comments:

Slide Evaluation

General Considerations

PD-L1 IHC 22C3 pharmDx evaluation should be performed by a qualified pathologist using a light microscope of diagnostic quality. Details of the PD-L1 IHC 22C3 pharmDx interpretation guidelines are reviewed on page 30. Before examining the patient specimen for PD-L1 staining, it is important to examine the controls to assess staining quality.

PD-L1 interpretation is best assessed by requesting 3 serial tissue sections (H&E, PD-L1 stain, and NCR stain) so that if the H&E is first assessed and is acceptable, IHC staining of the remaining 2 serial sections is likely to be acceptable.

Each PD-L1 IHC 22C3 pharmDx is configured with Control Cell Line Slides that should be included in each IHC run. Guidelines on interpreting the Control Cell Line Slide are reviewed to the right. In-house control tissue slides should also be assessed with every IHC run.

Specimen Adequacy

Confirm the Presence of at Least 100 Viable Tumor Cells

A hematoxylin and eosin (H&E) stained section is recommended for the evaluation of specimen adequacy. PD-L1 IHC 22C3 pharmDx and the H&E staining should be performed on serial sections from the same paraffin block of the specimen.

A minimum of 100 viable tumor cells must be present in the PD-L1 stained slide for the specimen to be considered adequate for PD-L1 evaluation.

Instructions for Patient Specimens With Less Than 100 Viable Tumor Cells

Tissue from a deeper level of the block, or potentially another block, could have a sufficient number of viable tumor cells for PD-L1 IHC 22C3 pharmDx testing.

Evaluating Controls




Figure 6: Each Control Cell Line Slide contains sections of cell pellets with positive and negative PD-L1 expression.

PD-L1 IHC 22C3 pharmDx Control Cell Line Slide

Examine the PD-L1 IHC 22C3 pharmDx Control Cell Line Slide to determine that reagents are functioning properly. Each slide contains sections of cell pellets with positive and negative PD-L1 expression (Figure 6). Assess the percentage of positive cells and the staining intensity. If any staining of the Control Cell Line Slide is not satisfactory, all results with the patient specimens should be considered invalid.

Evaluate the overall staining intensity using the following guide:

0	Negative
1+	Weak intensity
2+	Moderate intensity
3+	Strong intensity

Positive Control Cell Pellet

The following staining is acceptable for the PD-L1 positive cell pellet (Figure 7):

- Cell membrane staining of $\geq 70\%$ of cells
- $\geq 2+$ average staining intensity
- Non-specific staining $< 1+$ intensity




Figure 7: Positive cell pellet with acceptable staining of PD-L1 IHC 22C3 pharmDx Control Cell Line Slide (20x magnification).

Negative Control Cell Pellet

For the PD-L1 negative cell pellet, the following staining is acceptable (Figure 8):

- The majority of cells should demonstrate no staining. Note: The presence of 10 or fewer cells with distinct cell membrane staining is acceptable
- Any background staining is less than 1+ staining intensity




Figure 8: Negative cell pellet with no staining of PD-L1 IHC 22C3 pharmDx Control Cell Line Slide (20x magnification).

Positive and Negative In-house Control Tissue (Gastric Cancer)

Examine the positive in-house gastric cancer control tissue to determine that the tissues are correctly prepared and reagents are functioning properly. The ideal positive control tissue provides a complete dynamic representation of weak to moderate staining of tumor cells and tumor-associated mononuclear inflammatory cells (MICs) (Figure 9). If staining of positive in-house control tissue is not satisfactory, all results with the patient specimen should be considered invalid.




Figure 9: Ideal positive in-house control tissue (10x magnification).

The ideal negative control tissue should demonstrate no staining on tumor cells and immune cells (Figure 10). However, because prevalence of PD-L1 expression on immune cells is high, positively staining immune cells are acceptable. Examine the negative in-house control tissue to determine the expected staining. The variety of different cell types present in most tissue sections offers internal negative control sites; this should be verified by the user.

If unwanted staining occurs in the in-house control tissues, results with the patient specimen should be considered invalid.




Figure 10: Ideal negative in-house control tissue demonstrating lack of staining (10x magnification).

Note: Gastric cancer specimens include gastric or gastroesophageal junction adenocarcinoma.

Optional Control Tissue

FFPE tonsil may also be used as an optional control specimen. Tonsil stained with PD-L1 should exhibit strong membrane staining in portions of the crypt epithelium and weak to moderate membrane staining of the follicular macrophages in the germinal centers (Figure 11).

PD-L1 expression of the endothelium, fibroblasts, and the surface epithelium should be absent.




Figure 11: Tonsil stained with PD-L1 primary antibody exhibiting strong membrane staining in portions of the crypt epithelium (A) and weak to moderate membrane staining of follicular macrophages in the germinal centers (B) (10x magnification).

Do not use in-house control tissue as an aid in interpretation of patient results.

Negative Control Reagent (NCR)

Examine the slides stained with the NCR to identify non-specific background staining that may interfere with PD-L1 staining interpretation, making the specimen non-evaluable. Satisfactory performance is indicated by the absence of staining (Figure 12).




Figure 12: Ideal negative in-house control tissue stained with Negative Control Reagent (10x magnification).

Negative Control Reagent stained slides indicate non-specific background staining and allow better interpretation of patient specimens stained with the primary antibody.

Slide Evaluation Flowchart




Figure 13: Recommended order of slide evaluation.

Combined Positive Score

Definition of Combined Positive Score (CPS)

PD-L1 expression in gastric or GEJ adenocarcinoma is determined by using Combined Positive Score (CPS), which is the number of PD-L1 staining cells (tumor cells, lymphocytes, macrophages*) divided by the total viable tumor cells, multiplied by 100. Although the result of the calculation can exceed 100, the maximum score is defined as CPS 100.

CPS is defined accordingly:

$$\text{CPS} = \frac{\text{\# PD-L1 staining cells (tumor cells, lymphocytes, macrophages)}}{\text{Total \# viable tumor cells}} \times 100$$

* Macrophages and histiocytes are considered the same cells

CPS Numerator Inclusion and Exclusion Criteria

Any convincing partial or complete linear membrane staining ($\geq 1+$) of viable tumor cells that is perceived as distinct from cytoplasmic staining is considered PD-L1 staining and should be included in scoring.

Any convincing membrane and/or cytoplasmic staining ($\geq 1+$) of lymphocytes and macrophages (mononuclear inflammatory cells, MICs) within tumor nests and/or adjacent supporting stroma is considered PD-L1 staining and should be included in scoring. Only MICs directly associated with the response against the tumor are scored.

See Table 1 on page 24 for additional CPS numerator inclusion/exclusion criteria.

Determining Combined Positive Score

- At low magnifications (4x, 10x), examine all well-preserved tumor areas. Evaluate overall areas of PD-L1 staining and non-staining tumor cells, keeping in mind that partial membrane staining or 1+ membrane staining may be difficult to see at low magnifications. Ensure there are at least 100 viable tumor cells in the sample
 - A minimum of 100 viable tumor cells must be present in the PD-L1 stained slide (biopsy and resection) for the specimen to be considered adequate for evaluation
- If patient specimens include more than one biopsy (ie. 3-5 endoscopic biopsies) on a slide, all tissue on the slide needs to be evaluated to generate a single CPS for determining the PD-L1 expression level. Each biopsy should not be reported independently
- For specimens with less than 100 viable tumor cells, tissue from a deeper level of the block or potentially another block could have a sufficient number of tumor cells for evaluation of PD-L1 expression

- At higher magnification (20x), evaluate PD-L1 expression and calculate CPS:
 - o Determine the total number of viable tumor cells, both PD-L1 staining and non-staining (CPS denominator)
 - o Determine the number of PD-L1 staining cells (tumor cells, lymphocytes, macrophages) (CPS numerator; see Table 1 for additional CPS numerator inclusion/exclusion criteria)
 - o Calculate CPS
- Evaluation of membrane staining at 20x magnification is recommended

Table 1: CPS numerator inclusion/exclusion criteria

Tissue Elements	Included in the Numerator	Excluded from the Numerator
Tumor Cells	Convincing partial or complete linear membrane staining (at any intensity) of viable invasive gastric or GEJ adenocarcinoma tumor cells	<ul style="list-style-type: none"> - Non-staining tumor cells - Tumor cells with only cytoplasmic staining - Adenoma, dysplasia, and carcinoma in situ
Immune Cells	Membrane and/or cytoplasmic* staining (at any intensity) of mononuclear inflammatory cells (MICs) within tumor nests and adjacent supporting stroma†: <ul style="list-style-type: none"> - Lymphocytes (including lymphocyte aggregates) - Macrophages Only MICs directly associated with the response to the tumor are scored	<ul style="list-style-type: none"> - Non-staining MICs - MICs associated with adenoma, dysplasia, and carcinoma in situ - MICs (including lymphoid aggregates) associated with ulcers, chronic gastritis, and other processes not associated with the tumor - MICs associated with normal structures - Neutrophils, eosinophils, and plasma cells
Other Cells	Not included	<ul style="list-style-type: none"> - Normal cells - Stromal cells (including fibroblasts) - Necrotic cells and/or cellular debris

* In MICs, membrane and cytoplasmic staining are often indistinguishable due to high nuclear to cytoplasmic ratio. Therefore, membrane and/or cytoplasmic staining of MICs is included in the CPS numerator

† Adjacent MICs are defined as being within the same 20x field as the tumor. However, MICs that are NOT directly associated with the response to the tumor should be excluded

‡ Macrophages and histiocytes are considered the same cells

The CPS denominator includes all viable invasive tumor cells (PD-L1 staining and non-staining). Dysplasia, carcinoma in situ, and all other cells are excluded.

Suggested Methods

Agilent recommends that scoring be performed within the context of the pathologist’s past experience and best judgment in interpreting IHC stains. We offer three different examples of techniques that may be used when determining the respective Combined Positive Scores (CPS) of various staining patterns.

The entire IHC slide should be reviewed to determine which of the following example techniques may be used.


Example 1: Calculation of Combined Positive Score in a Small Tumor Area With Staining

At lower magnifications (4x, 10x): Evaluate the tumor area for convincing staining as described in “Determining Combined Positive Score” on page 23.

Assessment: 10% of area with staining, 90% of area without staining

At higher magnification (20x): Confirm there is no staining in areas that appeared void of staining at lower magnifications. Evaluate the area of staining to estimate the number of PD-L1 staining cells (tumor cells, lymphocytes, macrophages). Also estimate the total number of viable tumor cells (PD-L1 staining and non-staining tumor cells).

Assessment: There are approximately 100 viable tumor cells and about 80 PD-L1 staining cells (per the CPS numerator)



Calculate the Combined Positive Score of the entire tumor area:

Assessment:

CPS of area with staining:

$$\text{CPS} = \frac{\# \text{ PD-L1 staining cells}^*}{\text{Total \# viable tumor cells}} \times 100 = \frac{\sim 80 \text{ PD-L1 staining cells}}{100 \text{ tumor cells}} \times 100 = 80$$

CPS of entire tumor area: 10% × 80 = ~CPS 8

- PD-L1 staining tumor cell
- PD-L1 non-staining tumor cell
- PD-L1 staining mononuclear inflammatory cell (MIC)
- PD-L1 non-staining mononuclear inflammatory cell (MIC)

Clinical Interpretation: PD-L1 expression

* Including tumor cells, lymphocytes, macrophages


Figure 14: Example of tumor with small staining area.

Example 2: Calculation of Combined Positive Score in a Heterogeneous Tumor Area

At lower magnifications (4x, 10x): Visually divide the tumor area into regions with equal numbers of tumor cells.

At higher magnification (20x): Observe each region and estimate the total number of viable tumor cells and PD-L1 staining cells (tumor cells, lymphocytes, macrophages). Calculate the Combined Positive Score for each region.

Assessment: The four sections have ~80, ~30, ~50, and ~100 PD-L1 staining cells (tumor cells, lymphocytes, macrophages). Each section has a total of 100 tumor cells (including PD-L1 staining cells). The CPS for each section: ~CPS 80, ~CPS 30, ~CPS 50, and ~CPS 100



Calculate the Combined Positive Score of the entire tumor area:

Assessment:

Combined Positive Score:
 $(80 + 30 + 50 + 100) / 4 = \sim\text{CPS } 65$

$$\text{CPS} = \frac{\# \text{ PD-L1 staining cells (tumor cells, lymphocytes, macrophages)}}{\text{Total \# viable tumor cells}} \times 100$$

Clinical Interpretation: PD-L1 expression


Figure 15: Example with heterogeneous tumor area.

Example 3: Calculation of Combined Positive Score for a Near Cut-off Specimen (CPS 0–10)

At lower magnifications (4x, 10x): Evaluate the specimen for convincing staining as described in “Determining Combined Positive Score” on page 23.

At higher magnification (20x): Confirm that there is no staining in areas that appeared void of staining at lower magnifications. Evaluate all staining areas and estimate the total number of PD-L1 staining cells (tumor cells, lymphocytes, macrophages). Then re-evaluate the entire specimen (staining and non-staining areas) and estimate the total number of viable tumor cells (PD-L1 staining and non-staining tumor cells). Calculate the Combined Positive Score.

Assessment: Four areas of the tumor specimen have convincing staining. There are 8 PD-L1 staining cells (tumor cells, lymphocytes, macrophages) in the four staining areas. There are approximately 200 viable tumor cells present in the entire specimen



Calculate the Combined Positive Score of the entire tumor area:

Assessment:

Combined Positive Score:

$$\text{CPS} = \frac{\# \text{ PD-L1 staining cells}^*}{\text{Total \# viable tumor cells}} \times 100 = \frac{8 \text{ PD-L1 staining cells}}{200 \text{ tumor cells}} \times 100 = \text{CPS } 4$$

Clinical Interpretation: PD-L1 expression



* Including tumor cells, lymphocytes, macrophages

Figure 16: Example of near cut-off specimen (CPS 0–10).

Interpretation of CPS

The Combined Positive Score describes the PD-L1 expression of the specimen. See the table below for scoring guideline examples.

Table 2: CPS and PD-L1 expression

CPS	Expression Level	Image (20×)
< 1	No PD-L1 Expression	
≥ 1	PD-L1 Expression	

Identifying Patients With Urothelial Carcinoma for Treatment

PD-L1 IHC 22C3 pharmDx is the only companion diagnostic indicated as an aid in identifying patients with gastric or GEJ adenocarcinoma for treatment with KEYTRUDA® (pembrolizumab).

Clinical Validation of PD-L1 IHC 22C3 pharmDx in Previously Treated* Patients with Gastric or GEJ Adenocarcinoma

The clinical validity of PD-L1 IHC 22C3 pharmDx in evaluating PD-L1 expression (CPS \geq 1) in previously treated patients* with gastric or GEJ adenocarcinoma is based on the KEYTRUDA KEYNOTE-059 study sponsored by Merck Sharp & Dohme Corp. Specimens from previously treated patients with gastric or GEJ adenocarcinoma were tested for PD-L1 expression using PD-L1 IHC 22C3 pharmDx. Fifty-eight percent of enrolled patients had tumors that expressed PD-L1 with a Combined Positive Score (CPS) of greater than or equal to 1 (CPS \geq 1). Clinical efficacy of KEYTRUDA treatment in patients is presented in the Clinical Performance Evaluation section on page 68.

* At least 2 prior systemic treatments for advanced disease

Table 3: PD-L1 Prevalence^a in Patients with Gastric or GEJ Adenocarcinoma^b Enrolled in KEYNOTE-059^c

PD-L1 Expression	CPS < 1	CPS \geq 1
Prevalence (n)	42.0% (109)	58.0% (148)

a. Merck & Co., data on file.

b. Patients enrolled in KEYNOTE-059 Cohort 1.

c. International phase 2 study of pembrolizumab in patients with recurrent or metastatic gastric or GEJ adenocarcinoma who have experienced disease progression after two lines of prior therapy. ClinicalTrials.gov number NCT02335411.

Table 4: Tumor PD-L1 Expression by Specimen Type

Tumor Tissue	PD-L1 Expression (CPS \geq 1), n (%)	No PD-L1 Expression (CPS < 1), n (%)
Overall Study, n=257 (%)	148 (58)	109 (42)
Archival Tissue [†] , n=167	82 (49)	85 (51)
Newly Obtained Tissue [†] , n=90	66 (73)	24 (27)

[†] In the context of clinical trial KN059, a newly obtained biopsy was defined as a specimen obtained up to 6 weeks (42 days) prior to initiation of treatment on Day 1 (Cycle 1) with KEYTRUDA and with no additional anti-cancer treatment having been given after the specimen was obtained. Specimens that were > 42 days were classified as archival.

NOTE: If PD-L1 expression is not detected in an archival gastric or GEJ adenocarcinoma specimen, evaluate the feasibility of obtaining an additional tumor biopsy for PD-L1 testing.

PD-L1 IHC 22C3 pharmDx Testing Scheme

Use the following flow chart to help you understand which patients are indicated for treatment with KEYTRUDA® based on their CPS and treatment history.




Figure 17: Testing algorithm for PD-L1 IHC 22C3 pharmDx.

Reporting Results

Suggested information to include when reporting results with PD-L1 IHC 22C3 pharmDx.

PD-L1 IHC 22C3 pharmDx Summary of Sample Tested

Date of Run: _____

PD-L1 IHC 22C3 pharmDx Lot: _____

Staining Run Log ID: _____

Specimen ID: _____

Patient Identifiers: _____

Type of Service: IHC Stain with Manual Interpretation

Other: _____

PD-L1 Included in Gastric or GEJ Adenocarcinoma Comprehensive Panel: Yes: No:

Type of Tissue: Gastric Adenocarcinoma: Gastroesophageal Junction Adenocarcinoma: Other:

Number of Biopsies in Tissue Block*: _____

*If patient specimens include more than one biopsy (ie. 3–5 endoscopic biopsies) on a slide, all tissues on the slide need to be evaluated to generate a single CPS for determining the PD-L1 expression level. Each biopsy should not be reported independently.

PD-L1 Testing Results

Control Cell Line Slide Results: Pass: Fail:

Adequate Tumor Cells Present (≥ 100 cells): _____

PD-L1 IHC 22C3 pharmDx Result to Treating Physician

Combined Positive Score: _____

CPS ≥ 1 (PD-L1 expression): CPS < 1 (No PD-L1 expression):

Comments to Treating Physician:

- KEYTRUDA (pembrolizumab) is indicated for the treatment of patients with recurrent locally advanced or metastatic gastric or gastroesophageal junction adenocarcinoma with disease progression on or after two or more prior lines of therapy whose tumors have PD-L1 expression [Combined Positive Score (CPS) greater than or equal to 1] as determined by an FDA-approved test. See KEYTRUDA prescribing information for details
- In the context of clinical trial KN059, a newly obtained biopsy was defined as a specimen obtained up to 6 weeks (42 days) prior to initiation of treatment on Day 1 (Cycle 1) with KEYTRUDA and with no additional anti-cancer treatment having been given after the specimen was obtained. Specimens that were > 42 days were classified as archival

NOTE: If PD-L1 expression is not detected in an archival gastric or GEJ adenocarcinoma specimen, evaluate the feasibility of obtaining an additional tumor biopsy for PD-L1 testing.

Combined Positive Score Summary and Examples

Key Considerations in Scoring PD-L1 IHC 22C3 pharmDx Stained Specimens

By definition, PD-L1 staining cells in gastric or GEJ adenocarcinoma are:

- Tumor cells with convincing partial or complete linear membrane staining (at any intensity) that is perceived distinct from cytoplasmic staining
- Lymphocytes and macrophages (mononuclear inflammatory cells, MICs) within the tumor nests and/or adjacent supporting stroma with convincing membrane and/or cytoplasmic staining (at any intensity). MICs must be directly associated with the response against the tumor

PD-L1 expression status in gastric or GEJ adenocarcinoma is determined by Combined Positive Score (CPS), which is the number of PD-L1 staining cells (tumor cells, lymphocytes, macrophages) divided by the total number of viable tumor cells, multiplied by 100.

$$\text{CPS} = \frac{\text{\# PD-L1 staining cells (tumor cells, lymphocytes, macrophages)}}{\text{Total \# viable tumor cells}} \times 100$$

This section will define and illustrate all scoring inclusions and exclusions for accurate determination of Combined Positive Score. All images are gastric or GEJ adenocarcinoma unless otherwise noted in figure caption.

Image Guide for Interpretation of PD-L1 IHC 22C3 pharmDx Staining in Gastric or GEJ Adenocarcinoma

PD-L1 Staining Cells Included in the CPS Numerator

Tumor cells, lymphocytes, and macrophages exhibiting appropriate PD-L1 expression are defined as PD-L1 staining cells. All PD-L1 staining cells are included in the CPS numerator for determination of the Combined Positive Score (see Table 1 on page 24 for additional CPS numerator inclusion/exclusion criteria). Below are common staining characteristics of PD-L1 staining cells that must be included in the CPS numerator. All images are gastric or GEJ adenocarcinoma unless otherwise noted in figure caption.

Linear Membrane Staining

Tumor cells and tumor-associated mononuclear inflammatory cells (MICs: lymphocytes, macrophages) exhibiting convincing partial and/or complete linear membrane staining are considered PD-L1 staining cells. Convincing linear membrane staining can be present at any intensity and must be convincing at a 20x magnification.




Figure 18a: PD-L1 primary antibody exhibiting linear membrane staining of tumor cells (20x magnification).




Figure 18b: PD-L1 primary antibody exhibiting linear membrane staining of tumor-associated mononuclear inflammatory cells (A: lymphocytes, B: macrophages) (20x magnification). Note: Tumor-associated MICs are present within the tumor nests and/or adjacent supporting stroma and are directly associated with the response against the tumor.

Key point

Convincing linear membrane staining of tumor cells and convincing linear membrane and/or cytoplasmic staining of tumor-associated lymphocytes and macrophages should be included in the CPS numerator

Partial and Complete Linear Membrane Staining

Tumor cells and tumor-associated mononuclear inflammatory cells (MICs: lymphocytes, macrophages) can exhibit partial and/or complete linear membrane staining. Any partial or complete linear membrane staining observed at any intensity and observed at a 20x magnification must be included in the CPS numerator.

Note: Tumor-associated MICs are present within the tumor nests and/or adjacent supporting stroma and are directly associated with the response against the tumor.




Figure 19: PD-L1 primary antibody exhibiting partial (A) and complete (B) linear membrane staining of tumor cells (20x magnification).

Key point

Convincing partial and/or complete staining of tumor cells and tumor-associated lymphocytes and macrophages should be included in the CPS numerator

Weak Linear Membrane Staining

Tumor cells and tumor-associated mononuclear inflammatory cells (MICs: lymphocytes, macrophages) must exhibit convincing staining at any intensity, including weak 1+ intensity, at a 20x magnification.

Note: Tumor-associated MICs are present within the tumor nests and/or adjacent supporting stroma and are directly associated with the response against the tumor.




Figure 20a: PD-L1 primary antibody exhibiting weak but perceptible and convincing membrane staining of tumor cells (arrows) (20x magnification).




Figure 20b: PD-L1 primary antibody exhibiting weak but perceptible and convincing membrane staining of tumor-associated mononuclear inflammatory cells (arrows) (20x magnification).

Key point

Weak 1+ convincing staining of tumor cells and tumor-associated lymphocytes and macrophages should be included in the CPS numerator

Linear Membrane and Cytoplasmic Staining: Tumor Cells

Tumor cells with both convincing linear membrane staining ($\geq 1+$ intensity) and cytoplasmic staining at 20x magnification should be included in the CPS numerator. Tumor cells exhibiting only cytoplasmic staining are excluded from the CPS numerator, as this is considered non-specific staining.




Figure 21: PD-L1 primary antibody exhibiting linear membrane staining distinct from cytoplasmic staining (arrows) (20x magnification).

Key point

Tumor cells exhibiting convincing linear membrane staining that is distinct from cytoplasmic staining is included in the CPS numerator

Membrane and Cytoplasmic Staining: Tumor-associated Mononuclear Inflammatory Cells (MICs)

Tumor-associated lymphocytes and macrophages (mononuclear inflammatory cells, MICs) exhibiting convincing membrane and/or cytoplasmic staining at a 20x magnification ($\geq 1+$ intensity) are considered PD-L1 staining cells and should be included in the CPS numerator. Tumor-associated MICs are present within the tumor nests and/or adjacent supporting stroma and are directly associated with the response against the tumor.

Note: PD-L1 staining lymphocytes often have indistinguishable membrane and cytoplasmic staining due to a high nuclear to cytoplasmic ratio; PD-L1 staining macrophages often have distinct membrane staining and low cytoplasmic staining. All PD-L1 staining tumor-associated MICs should be included in the CPS numerator.




Figure 22a: PD-L1 primary antibody exhibiting linear membrane and/or cytoplasmic staining of tumor-associated mononuclear inflammatory cells (20x magnification).




Figure 22b: PD-L1 primary antibody exhibiting linear membrane and/or cytoplasmic staining of tumor-associated lymphocytes (20x magnification).




Figure 22c: PD-L1 primary antibody exhibiting linear membrane staining of tumor-associated macrophages (arrows) (20x magnification).

Key point

Tumor-associated lymphocytes and macrophages with convincing membrane and/or cytoplasmic staining should be included in the CPS numerator

Lattice and Interface Staining Patterns

PD-L1 staining tumor cells and tumor-associated mononuclear inflammatory cells can exhibit two distinct staining patterns: lattice and interface. Gastric or GEJ adenocarcinoma specimens may exhibit one or both patterns.

The lattice pattern of PD-L1 staining cells may be partial or complete, and present within tumor nests and/or adjacent supporting stroma. Staining can be present at one or several intensities.

Specimens with PD-L1 staining at the leading edge and/or margin of the tumor nests exhibit the interface pattern. PD-L1 staining can also be present in the adjacent supporting stroma. PD-L1 staining can be present at one or several intensities.




Figure 23a: PD-L1 primary antibody exhibiting the lattice pattern (20x magnification).




Figure 23b: PD-L1 primary antibody exhibiting the interface pattern (arrows) (10x magnification).

Heterogeneous Staining Intensities

Convincing staining of tumor cells (linear membrane) and tumor-associated lymphocytes and macrophages (membrane and/or cytoplasmic) is often heterogeneous, with various staining intensities present. At 20x magnification, any convincing staining of tumor cells and tumor-associated lymphocytes and macrophages at any intensity should be included in the CPS numerator.




Figure 24a: PD-L1 primary antibody exhibiting heterogeneous staining intensities of tumor cells (A: 1+ intensity, B: 2+ intensity, C: 3+ intensity) (20x magnification).




Figure 24b: PD-L1 primary antibody exhibiting heterogeneous staining intensities of tumor-associated mononuclear inflammatory cells (A: 1+ intensity, B: 2+ intensity, C: 3+ intensity) (20x magnification).

Key point

Convincing staining of tumor cells and tumor-associated lymphocytes and macrophages at all intensities should be included in the CPS numerator

Granular Staining

Tumor cells can exhibit a granular membrane staining pattern where membrane and cytoplasmic staining are indistinguishable. Only convincing staining of tumor cells ($\geq 1+$ intensity) observed at a 20x magnification should be included in the CPS numerator.




Figure 25: PD-L1 primary antibody exhibiting a granular membrane staining pattern (arrows) (20x magnification).

Key point

Granular staining of tumor cells must exhibit a convincing linear membrane pattern to be included in the CPS numerator

Diffuse Signet Ring Cell Carcinoma

Signet ring cells are commonly present in diffuse-type gastric or GEJ adenocarcinoma. The tumor cell nuclei are small, crescent-shaped, and present at the periphery of the cells. Signet ring cells exhibiting any convincing partial and/or complete membrane staining at a 20x magnification ($\geq 1+$ intensity) are considered PD-L1 staining cells and should be included in the CPS numerator. Signet ring cells exhibiting only cytoplasmic staining should be considered non-staining, as this is considered non-specific staining.

Key point

Convincing partial and/or complete membrane staining of signet ring cells should be included in the CPS numerator

Cells Excluded from CPS

Only tumor cells exhibiting PD-L1 membrane staining and MICs exhibiting PD-L1 membrane and/or cytoplasmic staining should be included in the CPS numerator. Below are other cells that can exhibit PD-L1 expression but should be excluded from the CPS calculation (CPS numerator or denominator).

Tumor Cells with Only Cytoplasmic Staining

Tumor cells exhibiting only cytoplasmic staining are not to be included in the CPS numerator, as this is considered non-specific staining. They should, however, still be included in the CPS denominator.




Figure 26: PD-L1 primary antibody exhibiting only cytoplasmic staining (arrows) (20x magnification).

Key point

Tumor cells exhibiting only cytoplasmic staining should not be included in the CPS numerator

Other Immune Cells Excluded from CPS

Various types of immune cells can exhibit PD-L1 staining, but only tumor-associated lymphocytes and macrophages should be included in the CPS calculation. Neutrophils, eosinophils, and plasma cells should be excluded from the CPS calculation.

Key point

Neutrophils, eosinophils, and plasma cells should be excluded from the CPS calculations

Non-tumor-associated MICs

Mononuclear inflammatory cells (MICs: lymphocytes and macrophages) commonly exhibit PD-L1 staining in gastric or GEJ adenocarcinoma specimens. Only PD-L1 staining MICs that are tumor-associated (present within tumor nests and/or adjacent supportive stroma; directly associated with the response against the tumor) should be included in the CPS calculation.

PD-L1 staining MICs that are not tumor-associated must be excluded from the CPS calculation. Examples of non-tumor-associated MICs include those associated with adenoma, dysplasia, carcinoma in situ, ulceration, chronic gastritis, normal cells, and other processes not associated with the tumor.




Figure 27a: PD-L1 staining of mononuclear inflammatory cells associated with normal gastric glands (10x magnification).




Figure 27b: PD-L1 staining mononuclear inflammatory cells contributing to gastritis (10x magnification).




Figure 27c: PD-L1 staining mononuclear inflammatory cells associated with no tumor cells (10x magnification).




Figure 27d: PD-L1 staining mononuclear inflammatory cells associated with dysplasia (arrows) (10x magnification).

Key point

Non-tumor-associated MICs should be excluded from the CPS calculation

Other Cell Types Excluded from CPS

Various other tissue elements can exhibit PD-L1 staining, but only PD-L1 staining tumor cells and tumor-associated lymphocytes and macrophages should be included in the CPS numerator. Normal cells (including ganglion cells), stromal cells (including fibroblasts), and endocrine cells should be excluded from the CPS calculation.




Figure 28a: PD-L1 staining of benign gastric glands (20x magnification).




Figure 28b: PD-L1 staining of fibroblasts (arrows) (20x magnification).




Figure 28c: PD-L1 staining of ganglion cells (arrows) (20x magnification).




Figure 28d: PD-L1 staining of endocrine cells (20x magnification).

Key point

PD-L1 staining normal cells, stromal cells, and endocrine cells should be excluded from the CPS calculation

CPS < 1 Case Examples

Case 1: CPS < 1




Figure 29a: 10x magnification.




Figure 29b: 20x magnification.

Figure 29a–29b: PD-L1 antibody exhibiting CPS < 1 at 20x magnification.

Case 2: CPS < 1




Figure 30a: 10x magnification.




Figure 30b: 20x magnification.

Figure 30a–30b: PD-L1 antibody exhibiting CPS < 1 at 20x magnification.

Case 3: CPS < 1




Figure 31a: 10x magnification.




Figure 31b: 20x magnification.

Figure 31a–31b: PD-L1 antibody exhibiting CPS < 1 with non-specific staining at 20x magnification.

Case 4: CPS < 1




Figure 32a: 10x magnification.




Figure 32b: 20x magnification.

Figure 32a–32b: PD-L1 antibody exhibiting CPS < 1 at 20x magnification.

Case 5: CPS < 1




Figure 33a: 10x magnification.




Figure 33b: 20x magnification.

Figure 33a–33b: PD-L1 antibody exhibiting CPS < 1 at 20x magnification.

Case 6: CPS < 1




Figure 34a: 10x magnification.




Figure 34b: 20x magnification.

Figure 34a–34b: PD-L1 antibody exhibiting CPS < 1 at 20x magnification.

CPS \geq 1 Case Examples

Case 7: CPS \geq 1




Figure 35a: 10x magnification.




Figure 35b: 20x magnification.

Figure 35a–35b: PD-L1 antibody exhibiting CPS \geq 1 at 20x magnification.

Case 8: CPS \geq 1




Figure 36a: 10x magnification.




Figure 36b: 20x magnification.

Figure 36a–36b: PD-L1 antibody exhibiting CPS \geq 1 at 20x magnification.

Case 9: CPS \geq 1




Figure 37a: 10x magnification.




Figure 37b: 20x magnification.

Figure 37a–37b: PD-L1 antibody exhibiting CPS \geq 1 at 20x magnification.

Case 10: CPS \geq 1




Figure 38a: 10x magnification.




Figure 38b: 20x magnification.

Figure 38a–38b: PD-L1 antibody exhibiting CPS \geq 1 at 20x magnification.

Case 11: CPS \geq 1




Figure 39a: 10x magnification.




Figure 39b: 20x magnification.

Figure 39a–39b: PD-L1 antibody exhibiting CPS \geq 1 at 20x magnification.

Case 12: CPS \geq 1




Figure 40a: 10x magnification.




Figure 40b: 20x magnification.

Figure 40a–40b: PD-L1 antibody exhibiting CPS \geq 1 at 20x magnification.

Near Cut-off Case Examples (CPS 0–10)

Challenging Case 1: Near Cut-off (CPS 0–10)




Figure 41: PD-L1 antibody exhibiting CPS 2 at 20x magnification.

Challenging Case 2: Near Cut-off (CPS 0–10)




Figure 42: PD-L1 antibody exhibiting CPS 2 at 20x magnification with weak PD-L1 staining.

Challenging Case 3: Near Cut-off (CPS 0–10)




Figure 43: PD-L1 antibody exhibiting CPS < 1 and non-specific staining at 20x magnification.

Challenging Case 4: Near Cut-off (CPS 0–10)




Figure 44: PD-L1 antibody exhibiting CPS 2 and non-specific staining at 20x magnification.

Challenging Case 5: Near Cut-off (CPS 0–10)




Figure 45: PD-L1 antibody exhibiting CPS < 1 at 20x magnification.

Challenging Case 6: Near Cut-off (CPS 0–10)




Figure 46: PD-L1 antibody exhibiting CPS 1 at 20x magnification.

Artifacts

The following pages provide examples of artifacts you may see when staining with PD-L1 IHC 22C3 pharmDx.

Non-specific Background Staining

Background staining is defined as diffuse, non-specific staining of a specimen. It is caused by several factors. These factors include, but are not limited to, pre-analytic fixation and processing of the specimen, incomplete removal of paraffin from sections, and incomplete rinsing of slides during staining.

The use of fixatives other than neutral buffered formalin may be a source of background staining. Background staining with PD-L1 IHC 22C3 pharmDx is rare.

Possible Causes of Background

- Improper drying of slides; ensure slides remain wet with buffer while loading onto Autostainer Link 48 and prior to initiating run
- Improper deparaffinization procedure
- Incomplete rinsing of reagents from slides

The non-specific background staining of the NCR-stained test specimen is useful in determining the level of background staining in the positive test specimen. All specimens must have $\leq 1+$ non-specific background staining.

Key point

All specimens must have $\leq 1+$ non-specific background staining




Figure 47: Negative Control Reagent (NCR) exhibiting acceptable non-specific background staining (20x magnification).

Edge Artifact

Commonly, edge artifacts are linked to the pre-analytic handling of the tissue.

- Inadequate processing of thick tissue samples may mimic edge artifact by rendering the central portion of the tissue suboptimally fixed relative to the peripheral areas. In these circumstances, the immunoreactivity based on the suboptimal central portion may be mistakenly interpreted as non-staining as optimal fixation is only present at the periphery
- Increased staining is observed around the periphery of the tissue specimen, known as the “edge artifact”
- Edge artifacts can be due to drying of the tissue specimen prior to fixation or during the staining procedure
- If the positive reaction is only at the edge of the tissue section (i.e., a few cell layers of staining at the periphery and ending abruptly with penetration into the centrally located tumor), scoring at the edge of the tissue specimen should be avoided

Key point

Scoring of the edge of a specimen should be avoided if staining is inconsistent with the rest of the specimen




Figure 48: PD-L1 primary antibody exhibiting edge artifact staining; edge staining should be excluded from the scoring (5x magnification).

Crush Artifact

Crush artifact is closely related to edge artifact. The compression of the tissues along the edges of the specimen can produce a linear staining that has to be interpreted as an artifact.

- Inadvertent crushing of the tissue occasionally occurs during sectioning, resulting in morphologically distorted cellular architecture
- When compared to surrounding cells, stronger staining may be observed in crushed cells. Crushed cells typically demonstrate condensed nuclei. Crushed cells should be avoided in scoring

Key point

Scoring of crush artifact should be avoided if staining is inconsistent with entire specimen




Figure 49: PD-L1 primary antibody exhibiting crush artifact (20x magnification).

Necrosis

Necrosis can be described as morphological changes indicative of cell death with undefined cellular detail. Necrosis is often present in gastric or GEJ adenocarcinoma specimens and should be excluded from scoring.

Key point

Scoring of necrotic areas should be excluded from the CPS calculation




Figure 50: PD-L1 primary antibody exhibiting staining of necrosis; necrosis should not be scored (20x magnification).

Poor Fixation

Standardization of fixation is very important when using PD-L1 IHC 22C3 pharmDx. Suboptimal fixation of tissues may give erroneous results.

Key point

Proper fixation is important for accurate diagnosis




Figure 51: PD-L1 primary antibody exhibiting poor tissue fixation (10x magnification).

Troubleshooting Guide

Troubleshooting Guidelines for PD-L1 IHC 22C3 pharmDx

For further troubleshooting help, contact your local Agilent representative.

Problem	Probable Cause	Suggested Action
No staining of slides	Programming error	Verify that the PD-L1 IHC 22C3 pharmDx program was selected for programming of slides
	Lack of reaction with DAB+ Substrate-Chromogen Solution	Verify that DAB+ Substrate-Chromogen Solution was prepared properly
	Sodium azide in wash buffer	Use only Dako Wash Buffer (Code K8007)
	Degradation of Control Cell Line Slide	Check kit expiration date and kit storage conditions on outside of package
Weak staining of specimen slides	Inappropriate fixation method used	Ensure that only neutral buffered formalin fixative and approved fixation methods are used
	Insufficient reagent volume applied	Check size of tissue section and reagent volume applied
	Inappropriate wash buffer used	Use only Dako Wash Buffer, Code K8007
Weak staining of specimen slides or of the positive cell line pellet on the Agilent-provided Control Cell Line Slide	Inadequate target retrieval	Verify that the 3-in-1 pre-treatment procedure was correctly performed
	Inappropriate wash buffer used	Use only Dako Wash Buffer, Code K8007
Excessive background staining of slides	Paraffin incompletely removed	Verify that the 3-in-1 pre-treatment procedure was correctly performed
	Slides dried while loading onto Autostainer Link 48	Ensure slides remain wet with buffer while loading and prior to initiating run
	Non-specific binding of reagents to tissue section	Check for proper fixation of the specimen and/or the presence of necrosis
Tissue detached from slides	Use of incorrect microscope slides	Use Dako FLEX IHC Microscope Slides, Code K8020, or Fisherbrand Superfrost Plus charged slides
Excessively strong specific staining	Inappropriate fixation method used	Ensure that only approved fixatives and fixation methods are used
	Inappropriate wash buffer used	Use only Dako Wash Buffer, Code K8007
Target Retrieval Solution is cloudy in appearance when heated	When heated, the Target Retrieval Solution turns cloudy in appearance	This is normal and does not influence staining

Note: If the problem cannot be attributed to any of the above causes, or if the suggested corrective action fails to resolve the problem, please call Agilent Technical Support for further assistance. Additional information on staining techniques and specimen preparation can be found in Dako Education Guide: Immunohistochemical Staining Methods (available from Agilent).

Clinical Performance Evaluation

The efficacy of KEYTRUDA® (pembrolizumab) was investigated in KEYNOTE-059 (KN059), a multicenter, non-randomized, open-label multi-cohort trial that enrolled 259 patients with gastric or gastroesophageal junction adenocarcinoma who progressed on at least 2 prior systemic treatments for advanced disease. Previous treatment must have included a fluoropyrimidine and platinum doublet. HER2/neu positive patients must have previously received treatment with approved HER2/neu targeted therapy. Patients with active autoimmune disease or a medical condition that required immunosuppression or with clinical evidence of ascites by physical exam were ineligible.

Patients received KEYTRUDA 200 mg every 3 weeks until unacceptable toxicity or disease progression that was symptomatic, rapidly progressive, required urgent intervention, occurred with a decline in performance status, or was confirmed at least 4 weeks later with repeat imaging. Patients without disease progression were treated for up to 24 months. Assessment of tumor status was performed every 6 to 9 weeks. The major efficacy outcome measures were ORR according to RECIST 1.1, as assessed by blinded independent central review, and duration of response.

PD-L1 expression level for 259 patient tumor biopsy or resection tissue (167 archival and 90 newly obtained; refer to definition in Table 5) was determined using PD-L1 IHC 22C3 pharmDx. PD-L1 expression level for 2 samples was not evaluable. Overall, fifty-eight percent (148/257) of the patients had tumors that expressed PD-L1 with a Combined Positive Score (CPS) ≥ 1 . Seventy-three percent (66/90) of patients whose tumors were newly obtained for PD-L1 testing and 49% (82/167) in patients whose archival tumors were tested expressed PD-L1 at CPS ≥ 1 (Table 5).

Of 148 patients with PD-L1 expression at CPS ≥ 1 , 143 were assessed to be either microsatellite stable (MSS) tumor status or had undetermined MSI or MMR status. The baseline characteristics of these 143 patients were: median age 64 years (47% age 65 or older); 77% male; 82% White, 11% Asian; and ECOG PS of 0 (43%) and 1 (57%). Eighty-five percent had M1 disease and 7% had M0 disease. Fifty-one percent had two and 49% had three or more prior lines of therapy in the recurrent or metastatic setting.

For the 143 patients that have PD-L1 expression (CPS ≥ 1), the ORR was 13.3% (95% CI: 8.2, 20.0); 1.4% had a complete response and 11.9% had a partial response. Among the 19 responding patients, the duration of response ranged from 2.8+ to 19.4+ months, with 11 patients (58%) having responses of 6 months or longer and 5 patients (26%) having responses of 12 months or longer.

Table 5: Tumor PD-L1 Expression by Specimen Type

Tumor Tissue	PD-L1 Expression (CPS ≥ 1), n (%)	No PD-L1 Expression (CPS < 1), n (%)
Overall Study, n=257 (%)	148 (58)	109 (42)
Archival Tissue*, n=167	82 (49)	85 (51)
Newly Obtained Tissue*, n=90	66 (73)	24 (27)

**In the context of clinical trial KN059, a newly obtained biopsy was defined as a specimen obtained up to 6 weeks (42 days) prior to initiation of treatment on Day 1 (Cycle 1) with KEYTRUDA and with no additional anti-cancer treatment having been given after the specimen was obtained. Specimens that were > 42 days were classified as archival.*

References

- Keytruda [prescribing information]. Kenilworth, NJ: Merck & Co.; 2017.
- PD-L1 IHC 22C3 pharmDx [package insert]. Carpinteria, CA: Dako, Agilent Pathology Solutions; 2017.
- Herbst RS, Soria JC, Kowanetz M, et al. Predictive correlates of response to the anti-PD-L1 antibody MPDL3280A in cancer patients. *Nature*. 2014;515(7528):563-567.
- Tumeh PC, Harview CL, Yearley JH, et al. PD-1 blockade induces responses by inhibiting adaptive immune resistance. *Nature*. 2014;515(7528):568-571.
- Hanahan D, Weinberg RA. Hallmarks of cancer: the next generation. *Cell*. 2011;144(5):646-674.
- Sharpe AH, Freeman GJ. The B7-CD28 superfamily. *Nat Rev Immunol*. 2002;2(2):116-126.
- Keir ME, Butte MJ, Freeman GJ, Sharpe AH. PD-1 and its ligands in tolerance and immunity. *Annu Rev Immunol*. 2008;26:677-704.
- Freeman GJ, Long AJ, Iwai Y, et al. Engagement of the PD-1 immunoinhibitory receptor by a novel B7 family member leads to negative regulation of lymphocyte activation. *J Exp Med*. 2000;192(7):1027-1034.
- Dong H, Strome SE, Salomao DR, et al. Tumor-associated B7-H1 promotes T-cell apoptosis: a potential mechanism of immune evasion. *Nat Med*. 2002;8(8):793-800.
- Brahmer JR, Tykodi SS, Chow LQ, et al. Safety and activity of anti-PD-L1 antibody in patients with advanced cancer. *N Engl J Med*. 2012;366(26):2455-2465.
- Topalian SL, Hodi FS, Brahmer JR, et al. Safety, activity, and immune correlates of anti-PD-1 antibody in cancer. *N Engl J Med*. 2012;366(26):2443-2454.

Learn more:

<https://www.agilent.com/chem/PDL122C3>

For countries outside of the United States, see the local KEYTRUDA product label for approved indications and expression cutoff values to guide therapy.

This information is subject to change without notice.

© Agilent Technologies, Inc. 2019
Published, July 17, 2019
29219 2019JUL17

