

Minimum System Requirements

In case of Windows 10

CPU: Intel(R) Xeon(R) E3-1225 V5 @ 3.30GHz
Memory: 8GB
HDD: 500 GB

In case of Windows 7

CPU: Intel(R) Xeon(R) E3-1225 @3.10GHz
Memory: 4GB
HDD: 250 GB

Tested System

In case of Windows 10

PC model: HP Z240 SFF
CPU: Intel(R) Xeon(R) E3-1225 V5 @ 3.30GHz
Memory: 8GB
HDD: 500 GB
Resolution: 1920 x 1080
Printer: HP LaserJet M506 (English / Chinese)
Canon LBP251 (Japanese)
Network: Peer to Peer connection to instrument
Corporate LAN connection via a switching HUB
ICP-MS: 7900 Series, 7800 Series, 7700 Series, 8900 Series or 8800 Series
OS: Windows 10 Pro (64bit) Version 1803, OS Build 17134
Windows 10 Pro (64bit) Version 1703, OS Build 15063
Bootp: Agilent Bootp service A.01.01
Excel: No Excel (Note: Lab QC doesn't work if Excel is NOT installed.)
@@@ @@@ Excel 2016 (32bit)
Other Software Modules:
.Net framework 4.7.2
VSTO:
Visual Studio 2010 Tools for Office Runtime(x64)
Visual C++ Runtime:
Visual C++2005 Redistributable
Visual C++2005 Redistributable (x64)
Visual C++2008 Redistributable (x64) 9.0.30729.17
Visual C++2008 Redistributable (x86) 9.0.30729.17
Visual C++2010 Redistributable (x64) 10.0.40219
Visual C++2010 Redistributable (x86) 10.0.40219
Visual C++2012 Redistributable (x64) 11.0.61030
Visual C++2012 Redistributable (x86) 11.0.61030
Visual C++2013 Redistributable (x64) 12.0.30501
Visual C++2013 Redistributable (x86) 12.0.30501
Visual C++2015 Redistributable (x64) 14.0.23026
Visual C++2015 Redistributable (x86) 14.0.23026

In case of Windows 7

PC model: HP Z240 SFF, Z230 SFF, Z220 SFF and Z210 SFF

CPU: Z240: Intel(R) Xeon(R) E3-1225 V5 @ 3.30GHz
 Z230: Intel(R) Xeon(R) E3-1225 V3 @ 3.20GHz
 Z220: Intel(R) Xeon(R) E3-1225 V2 @ 3.20GHz
 Z210: Intel(R) Xeon(R) E3-1225 @ 3.10GHz
 Memory: 4GB
 HDD: Z240, Z230: 500 GB
 Z220, Z210: 250 GB
 Resolution: 1920 x 1050
 Printer: HP LaserJet M506 (English / Chinese)
 Canon LBP251 (Japanese)
 Network: Peer to Peer connection to instrument
 Corporate LAN connection via a switching HUB
 ICP-MS: 7900 Series, 7800 Series, 7700 Series, 8900 Series or 8800 Series
 OS: Windows 7 Professional SP1 (64bit) or Windows 7 Enterprise SP1
 (64bit)
 Bootp: Agilent Bootp service A.01.01
 Excel: No Excel (Note: Lab QC doesn't work if Excel is NOT installed.)
 @@@@ Excel 2013 (32bit)

Other Software Modules:

.Net framework 4.7.2
 VSTO:
 Visual Studio 2010 Tools for Office Runtime(x64)
 Visual C++ Runtime:
 Visual C++2005 Redistributable
 Visual C++2005 Redistributable (x64)
 Visual C++2008 Redistributable (x64) 9.0.30729.17 (Z230, Z220,
 Z210)
 Visual C++2008 Redistributable (x86) 9.0.30729.4148 (Z240)
 9.0.30729.17 (Z230, Z220,
 Z210)
 Visual C++2010 Redistributable (x64) 10.0.40219
 Visual C++2010 Redistributable (x86) 10.0.40219
 Visual C++2012 Redistributable (x64) 11.0.61030
 Visual C++2012 Redistributable (x86) 11.0.61030
 Visual C++2013 Redistributable (x64) 12.0.30501
 Visual C++2013 Redistributable (x86) 12.0.30501
 Visual C++2015 Redistributable (x64) 14.0.23026
 Visual C++2015 Redistributable (x86) 14.0.23026